## Sensibilissimo Test Droga Salivare

(7 Gruppi droga: oltre 25 sostanze)

CARATTERISTICHE

### CODICI PARAF:

912112265 ( Oratect III HM11 ) linea Blu x ME-AMP-COC-OPI-TH-PCP

912112253 ( Oratect III HM12 ) linea Blu x ME-AMP-COC-OPI-TH-BZ

912154150 ( Oratect III HMA11 + **Alcool** ) linea Blu x ME-AMP-COC-OPI-TH-PCP - alcool

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Acetaminophen	Chloroquine	Ethyl Morphine (except OP assay)	MDEA (except AM,ME assays)	Phenylephrine (except ME assay)
Acetoacetic acid lithium salt	d-Chlorpheniramine	Flunitrazepam (except BZ assay)	MDMA (except ME assay)	-Phenylethylamine (except AM assay)
Acetone	Chlorpromazine	Flurazepam (except BZ assay)	Methylphenidate	d,l-Phenylpropanolamine hydrochloride
6-Acetylcodeine (except OP assay)	Cholesterol	Furosemide	Morphine (except OP assay)	Prazepam (except BZ assay)
6-Acetylmorphine (except OP assay)	Clomipramine	Gentisic acid	Morphine-3-beta -D-Glucuronide (except OP assay)	Prednisolone
Acetylsalicylic acid	Clonazepam (except BZ assay)	Glucose	Nalidixic acid	Procaine
Albumin	Cocaine (except CO assay)	Glutethimide	Nalorphine (except OP assay)	Promazine
Allobarbital	Codeine (except OP assay)	Guaiacol Glyceryl Ether	Naloxone	Promethazine
Alphenal	Cortisone	Hemoglobin	Naltrexone	d-Propoxyphene
Alprazolam (except BZ assay)	I-Cotinine	Heroin (except OP assay)	hydrochloride	Protriptyline
Amitriptyline	Creatine	Hippuric acid	d-Naproxen	d-Pseudoephedrine
Amobarbital	Creatinine	Hydrochlorothizide	Niacinamide	HCl (except ME assay)
Amoxapine	Cyclobenzaprine	Hydrocodone (except OP assay)	Nitrazepam (except BZ assay)	Quinidine
Amoxicillin	Delorazepam (except BZ assay)	Hydrocortisone	11-Nor-Delta 8-THC-9-COOH (except TH assay)	Ranitidine
d-Amphetamine (except AM assay)	Deoxycortisone acetate	Hydromorphone (except OP assay)	11-Nor-Delta 9-THC-9-COOH (except TH assay)	Riboflavin
d,l-Amphetamine (except AM assay)	Desipramine (except ME assay)	11-Hydroxy-.A-9-Tetrahydrocannabinol (except TH assay)	Nordiazepam (except BZ assay)	Salicylic acid
l-Amphetamine (except AM assay)	Dextromethorphan	p-Hydroxymethamphetamine (Pholderin) (except ME assay)	Nordoxepin hydrochloride	Secobarbital
Ampicillin	Diazepam (except BZ assay)	3-Hydroxytyramine	d,l-Norephedrine hydrochloride	Serotonin
Apomorphine	Dihydrocodeine (except OP assay)	Ibuprofen	Norethindrone	Sodium Chloride
Aprobarbital	4Dimethylaminoantipyrine	Imipramine	d-Norpropoxyphene	Sulfamethazine
l-Ascorbic Acid	Diphenhydramine	d,l-Isoproterenol	Nortriptyline hydrochloride	Sulindac
Aspartame	Dopamine	l-Isoproterenol HCl	Oxalic Acid	Temazepam (except BZ assay)
Atropine	Doxepin hydrochloride	Lidocaine	Oxazepam (except BZ assay)	Tetracycline
Barbital	Doxylamine	Lorazepam (except BZ assay)	Oxolinic acid	.-8-Tetrahydrocannabinol (except TH assay)
Benzillic acid	Ecgonine (except CO assay)	Lormetazepam (except BZ assay)	Oxycodone	.-9-Tetrahydrocannabinol (except TH assay)
Benzocaine	Ecgonine Methyl Ester	Meperidine	Papaverine	Thiamine
Benzoic acid	d,l-Ephedrine (except ME assay)	d,l-Methadone	Penicillin-G (Benzylpenicillin)	Thioridazine
Benzoyllecgonine hydrate (except CO assay)	l-Ephedrine	d-Methamphetamine (except ME assay)	Pentazocine	Triazolam (except BZ assay)
Bilirubin	1R, 2S l- Ephedrine (except ME assay)	d,l-Methamphetamine (except ME assay)	Perphenazine	
Bromazepam (except BZ assay)	1S, 2R d-Ephedrine	l-Methamphetamine (except ME assay)	Phencyclidine (except PC assay)	
d-Brompheniramine	l-Epinephrine	Methaqualone	Pheniramine	
Buprenorphine	Erythromycin	Methoxyphenamine(except ME assay)	Phenobarbital	
Butalbital	Estazolam (except BZ assay)	2-Methylamine-Propiophenone HCl	Phenothiazine	
Butethal	β-Estradiol	MDA (except AM assay)	Phentermine (except AM assay)	
Caffeine	Estrone-3-sulfate			
Cannabidiol	potassium salt			
Cannabinol (except TH assay)	Ethanol			
Chloral Hydrate	Ethylidene-1,5-Dimethyl-1-3,3-Diphenylpyrrolidine			
Chlordiazepoxide (except BZ assay)	Perchlorate salt			
Chloroamphetamine (DL-p-) Hydrochloride (except AM assay)				

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



Drug Test	Approximate Concentration (ng/ml)	Approximate % Cross Reactivity
<b>ME/MDMA</b>		
Desipramine	10000	0.25%
d,l-Ephedrine	1000	2.5%
1R, 2S l-Ephedrine	1000	2.5%
p-Hydroxymethamphetamine	1000	2.5%
MDEA	300	8.3%
MDMA	25	100%
d,l-Methamphetamine	30	83%
d-Methamphetamine	25	100%
l-Methamphetamine	500	5%
Methoxyphenamine	2500	1%
Phenylephrine	5000	0.5%
d-Pseudoephedrine HCl	5000	0.5%
Trimethobenzamide	4000	0.6%
<b>TH</b>		
Cannabinol	80	50%
Δ-8-tetrahydrocannabinol	100	40%
Δ-9-tetrahydrocannabinol	40	100%
11-nor-Δ-8-THC-9-COOH	10	400%
11-nor-Δ-9-THC-9-COOH	10	400%
11-hydroxy-Δ9-THC	400	10%
<b>CO</b>		
Benzoylcegonine	18	110%
Cocaine	20	100%
Ecgonine	5000	0.4%
<b>AM</b>		
d-Amphetamine	25	100%
d,l-Amphetamine	40	62.5%
l-Amphetamine	800	3.2%
d,l-p-Chloramphetamine	200	12.5%
MDA	40	62.5%
MDEA	100	25%
Phentermine	100	25%
β-Phenylethylamine	8000	0.3%
Tyramine	8000	0.3%
<b>OP</b>		
6-Acetylcodeine	20	50%
6-Acetylmorphine	12	83%
Codeine	10	100%
Dihydrocodeine	10	100%
Ethyl morphine	60	17%
Heroin	15	67%
Hydrocodone	60	17%
Hydromorphone	70	14%
Morphine	10	100%
Morphine-3-beta-D-Glucuronide	25	40%
Nalorphine	100	10%
<b>PC</b>		
Phencyclidine	4	100%
<b>BZ</b>		
Alorazolam	4	125%
Bromazepam	4	125%
Chlordiazepoxide	50	10%
Clobazam	10	50%
Clonazepam	20	25%
Delorazepam	5	100%
Diazepam	5	100%
Estazolam	3	167%
Flunitrazepam	8	83%
Flurazepam	10	50%
Lorazepam	10	50%
Lormetazepam	15	33%
Nitrazepam	4	25%
Nordiazepam	3	67%
Oxazepam	5	100%
Frazeepam	10	50%
Temazepam	5	100%
Triazolam	10	50%

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